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1. McMichael, AJ ed: Leukocyte Typing III, White Cell Differentiation Antigens: 1987. Oxford: Oxford University Press, p 202, 206.
2. Reinherz, EL and Schlossman, SF: 1980. The differentiation and function of human T lymphocytes. Cell 19: 821-827.
3. Aiuta F, Cerottini J-C, Coombs RRA, Cooper M, Dickler HB, Froland S, Fudenberg HH, Greaves MF, Grey HM, Kunkel HG, Natvig J, Preud'homme J-L, Rabellino E, Ritts RE, Rowe DS, Seligmann M, Siegal FP, Stjernsward J, Terry WD and Wybran J: 1975. Identification, enumeration and isolation of B and T lymphocytes from human peripheral blood. International Union of Immunological Societies (IUIS), Report - July 1974. Clin Immunol and Immunopathol 3: 584-597
4. Reinherz EL, Haynes BF, Nadler LM and Bernstein ID: 1986. Leukocyte Typing II. New York: Springer-Verlag. Vol.2, p.8, 15-20, 37.
5. Reinherz, EL, Meuer, SC and Schlossman, SF: 1983. The delineation of antigen receptors on human T lymphocytes. Immunol Today 4: 5-8.
6. Reinherz, EL, Morimoto, C, Fitzgerald, KA, Hussey, RE, Daley, JF and Schlossman, SF: 1982. Heterogeneity of human T4+ inducer T cells defined by a monoclonal antibody that delineates two functional subpopulations. J Immunol 128: 463-468.
7. de Martini, RM and Parker, JW: 1989. Immunologic alterations in Human immunodeficiency virus infection: A review. J Clin Lab Anal 3: 56-70.
8. Morimoto, C, Letvin, NL, Distaso, JA, Aldrich, WR and Schlossman, SF: 1985. The isoform and characterization of the human suppressor inducer T cell subset. J Immunol 134: 1508-1515.
9. Morimoto, C, Letvin, NL, Distaso, JA, Brown, HM and Schlossman, SF: 1986. The cellular basis for the induction of antigen-specific T8-suppressor cells. Eur J Immunol 16: 198-204.
10. Meuer, SC, Schlossman, SF and Reinherz, EL: 1982. Clonal analysis of human cytotoxic T lymphocytes: T4+ and T8+ effector T cells recognize products of different major histocompatibility regions. Proc Natl Acad Sci USA 79: 4395-4399.
11. Reinherz, EL, Cooper, MD and Schlossman, SF: 1981. Abnormalities of T cell maturation and regulation in human beings with immunodeficiency disorders. J Clin Invest 68: 699-705.
12. Schmidt, RE: 1989. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut 59:200-206.
13. Collier, AC, Meyers, JD, Corey, L, Murphy, VL, Roberts, PL and Handsfield, HH: 1987. Cytomegalovirus infection in homosexual men. Relationship to sexual

- practices, antibody to human immunodeficiency virus, and cell-mediated immunity. Am J Med 82: 593-601.
14. Taylor, MGJ, Fahey, JL, Detels, R and Giorgi, JV: 1989. CD4 percentage, CD4 number and CD4:CD8 ratio in HIV infection: How to choose and how to use. J AIDS 2: 114-124.
  15. Morimoto, C, Hafler, DA, Weiner, HL, Letvin, NL, Hagan, M, Daley, J and Schlossman, SF: 1987. Selective loss of the suppressor-inducer T-cell subset in progressive multiple sclerosis. Analysis with anti-2H4 monoclonal antibody. N Engl J Med 316:67-72.
  16. Reinherz, EL, Weiner, HL, Hauser, SL, Cohen, JA, Distaso, JA and Schlossman, SF: 1980. Loss of suppressor T cells in active multiple sclerosis. Analysis with monoclonal antibodies. N Engl J Med 303: 125-129.
  17. Weiner, JL, Hafler, DA, Fallis, RJ, Johnson, D, Ault, KA and Hauser, SL: 1984. Altered blood T-cell subsets in patients with multiple sclerosis. J Neuroimmunol 6: 115-121.
  18. Smolen, JS, Chused, TM, Leiserson, WM, Reeves, JP, Alling, DW and Steinberg, AD: 1982. Heterogeneity of immunoregulatory T cell subsets in systemic lupus erythematosus. Am J Med 72:783-790.
  19. Smolen, JS, Morimoto, C, Stetnberg, AD, Wolf, A, Schlossman SF, Steinberg, RT, Penne, E, Reinherz, EL, Rerchlin, M and Chused, TM: 1985. Systemic lupus erythematosus: delineation of subpopulations by clinical, serologic and T cell marker analysis. Am J Med Sci 289: 139-147.
  20. Morimoto, C Reinherz, EL, Schlossman, SF, Shur, PH, Mills, JA and Steinberg, SD: 1980. Alterations in immunoregulatory T cell subsets in active systemic lupus erythematosus. J Clin Invest 66: 1171-1174.
  21. Raziuddin, S, Nur, MA and Alwabel, AA: 1989. Selective loss of the CD4+ inducers of suppressor T cell subsets (2H4+) in active systemic lupus erythematosus. J Rheumatol 16: 1315-1319.
  22. Sato, K, Miyasaka, N, Yamaoka, K, Okuda, M, Yata, J and Nishioka, K: 1987. Quantitative detect of CD4+2H4+ cells in systemic lupus erythematosus and Sjogren's syndrome. Arthritis Rheum 30: 1407-1411.
  23. Ramos, EL, Turka, LA, Leggat, JE, Wood, IG, Milford, EL and Carpenter, CB: 1989. Decrease in phenotypically defined T helper inducer cells (T4+4B4+) and increase in T suppressor effector cells (T8+2H4+) in stable renal allograft recipients. Transplantation 47: 465-471.
  24. Pedrazzini, AS, Freedman, AS, Andersen, J, Heftin, L, Anderson, K, Takvorian, R, Canellos, GP, Whitman, J, Coral, F, Ritz, J and Nadler, LM: 1989. Anti-B cell monoclonal antibody purged autologous bone marrow transplantation for B-cell non-Hodgkin's lymphoma: Phenotypic reconstitution and B-cell function. Blood 74: 2203-2211.
  25. Reinherz, EL, Kung, PC, Goldstein, G and Schlossman, SF: 1979. A monoclonal antibody with selective reactivity with functionally mature human thymocytes and all peripheral human T cells. J Immunol 123: 1312-1317.
  26. Reinherz, EL, Haynes, BF, Nadler, LM and Bernstein, ID, eds: Leukocyte, Typing II, Human T Lymphocytes: 1986. New York: Springer-Verlag, p 8-9.
  27. Koepke JA and Landay AL. Precision and accuracy of absolute lymphocyte counts. Clin Immunol Immunopathol, 1989; 52:19-27.

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