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- Schlossman SF, et al., eds: 1995. Leukocyte a Typing V. Oxford University Press. pp.262,263,268,270,507-511, 1398-1400.
- Bernard A, et al., eds: 1984. Leukocyte Typing. New York: Springer-Verlag. pp.28, 41-42,44,196.
- McMichael AJ, ed: 1987. Leukocyte a Typing III. Oxford: Oxford University Press. pp.38,40,42,43,116,167,170-172,176,199, 202,206,302-308,315,475,796-801.
- Reinherz EL and Schlossman SF: 1980. The differentiation and function of human T lymphocytes. Cell 19:821-827.
- Robertson MJ, et al.: 1996. Characterization of a cell line, NKL, derived from an aggressive human natural killer cell leukemia. Experimental Hematology 24:406-415.
- Aiuta F, et al.: 1975. Identification, enumeration and isolation of B and T lymphocytes from human peripheral blood, International Union of Immunological Societies (IUIS), Report - July 1974. Clin Immunol and Immunopathol 3: 584-597
- Foon KA and Todd RF: 1986. Immunologic classification of leukemia and lymphoma. Blood 68: 1-31.
- Hercend T, et al. 1982. Generation of a cloned NK cell line derived from the "null cell" fraction of human peripheral blood. J Immunol 129: 1299-1305.
- Drexler HG, et al.: 1988. Routine immunophenotyping of acute leukaemias. Blut 57:327-339.
- Robertson MJ and Ritz J: 1990. Biology and clinical relevance of human natural killer cells. Blood 76: 2421-2438.
- Reinherz EL, et al.: 1986. Leukocyte Typing II. New York: Springer-Verlag. V012, p.8, 15-25, 37.
- Newman W, et al.: 1984. Immunological and immunochemical aspects of the T-200 family of glycoproteins. Mol Immunol 21:1113-1121.
- Fabre JW, et al.: 1977. Quantitative serological analysis of a rabbit anti-rat serum and preliminary biochemical characterization of the major antigen recognized. Transplantation 23: 349-359.
- Coffman RL, et al.: 1981. B220: A B cell-specific member of the T200 glycoprotein family. Nature 289:681-683.
- Dalchau R, et al.: 1980. Identification with a monoclonal antibody of a predominantly B lymphocyte-specific determinant of the human leukocyte common antigen. J Exp Med 153:753-765.
- Omary MB, et al.: 1980. Human homologue of murine T200 glycoprotein. J Exp Med 152: 842-852.
- Dalchau R, et al.: 1981. Monoclonal antibody to a human leukocyte-specific membrane glycoprotein probably homologous to the leukocyte-common (L-C) antigen of the rat. Eur J Immunol 10: 737-744.
- Nadler LM, et al.: 1983. B4, human B lymphocyte associated antigen expressed on normal, mitogen activated and malignant B lymphocytes. J Immunol 131:244-250.
- Reinherz EL, et al.: 1983\* The delineation of antigen receptors on human T lymphocytes. Immunol Today 4: 5-8.
- Reinherz EL, et al.: 1982. Heterogeneity of human T4+ inducer T cells defined by a monoclonal antibody that delineates two functional subpopulations. J Immunol 128:463-468.
- de Martini RM and Parker JW: 1989. Immunologic alterations in human immunodeficiency virus infection: A review. J Clin Lab Anal 3: 56-70.
- Reinherz EL, et al.: 1981. Antibody directed at a surface structure inhibits cytotoxicity but not suppressor function of human T lymphocytes. Nature 294: 168-170.
- Caligiuri M, et al.: 1987. Phenotypic and functional deficiency of natural killer cells in patients with chronic fatigue syndrome. J Immunol 139: 3306-3313.
- Morimoto C, et al.: 1985. The isolation and characterization of the human suppressor inducer T cell subset. J Immunol 134: 1508-1515.

- Morimoto C, et al.: 1986. The cellular basis for the induction of antigen-specific T8-suppressor cells. Eur J Immunol 16:196-204.
- Meuer SC, et al.: 1982. Clonal analysis of human cytotoxic T lymphocytes: T4+ and T8+ effector T cells recognize products of different major histocompatibility regions. Proc Natl Acad Sci USA 79: 4395-4399.
- Barclay AN, et al., ed: 1993. The Leucocyte Antigen Facts Book. Academic Press, pp. 106-109.
- Tedder TF and Isaacs CM: 1989. isolation of a CDNAS encoding the CD19 antigen of human and mouse B lymphocytes: a new member of the immunoglobulin superfamily. J Immunol 143:712-717.
- Tedder TF, et al.: 1997. The CD19-CD21 complex regulates signal transduction thresholds governing humoral immunity and autoimmunity. Immunity 6: 107-118.
- Griffin JD, et al.: 1983. Characterization of an antigen expressed by human natural killer cells. J Immunol 130: 2947-2951.
- Hercend T, et al.: 1985. Generation of monoclonal antibodies to a human natural killer clone. Characterization of two natural killer-associated antigens, NKH1A and NKH2. J Clin Invest 75: 932-943.
- Shmidt RE, et al.: 1986. A subset of natural killer cells in peripheral blood displays a mature T cell phenotype. J Exp Med 164: 351-356.
- Schmidt RE, et al.: 1987. Enhancement of natural killer function through activation of the T11E rosette receptor. J Clin Invest 79:305-308.
- Reinherz EL, et al.: 1982. Antigen recognition by human T lymphocytes is linked to surface expression of the T3 molecular complex. Cell 30: 735-743.
- Meuer SC, et al.: 1983. Evidence for the T3-associated 90K heterodimer as the T-cell antigen receptor. Nature 303: 808-810.
- Benjamini E and Leskowitz S: 1991. Immunology: A Short Course. Second Edition. New York: Wiley-Liss. p.211-244.
- Reinherz EL, et al.: 1980. The cellular basis for viral-induced immunodeficiency: Analysis by monoclonal antibodies. J Immunol 125: 1269-1274.
- Felsenstein D, et al.: 1985. Phenotypic properties of atypical lymphocytes in cytomegalovirus-induced mononucleosis. J Infect Dis 152: 198-203.
- Rinaldo CR, et al.: 1983. Lymphocyte subsets and natural killer cell responses during cytomegalovirus mononucleosis. Infect Immun 40: 472-477.
- Goldstein G, et al.: 1982. Immunoregulatory changes in human disease detected by monoclonal antibodies to T lymphocytes. In: Monoclonal Antibodies in Clinical Medicine. McMichael AJ and Febre JW, eds. New York, NY: Academic Press. p.39-70.
- Schmidt RE: 1989. Monoclonal antibodies for diagnosis of immunodeficiencies. Bull 59: 200-206.
- Pedrazzini A, et al.: 1989. Anti-B cell monoclonal antibody purged autologous bone marrow transplantation for B-cell non-Hodgkin's lymphoma: Phenotypic reconstitution and B-cell function. Blood 74: 2203-2211.
- Preijers FWMB, et al.: 1989. Autologous transplantation of bone marrow purified in vitro with anti-CD7-(WT1)- Ricin A Immunotoxin in T-cell lymphoblastic leukemia and lymphoma. Blood 74: 1152-1158.
- Ramos EL, et al.: 1989. Decrease in phenotypically defined T helper inducer cells (T4+4B4+) and increase in T suppressor effector cells (T8+2H4+) in stable renal allograft recipients. Transplantation 47:465-471.
- Trinchieri G: 1989. Biology of natural killer cells. Adv Immunol 47:187.
- Fauci AS: 1988. The human deficiency virus: Infectivity and mechanism of pathogenesis. Science 239:617-622.
- Taylor MGJ, et al.: 1989. CD4 percentage. CD4 number and C04:CD8 ratio in HIV infection: How to choose and how to use. J AIDS 2:114-124.
- Center for Disease Control and Prevention: 1997. Revised guidelines for performing CD4+ T-cell determinations in persons infected with human immunodeficiency virus (HIV). MMWR 1997, 46 (No.RR.2): 1-29.
- Knapp W, et al., ed: 1989. Leukocyte a Typing IV. Oxford University Press. pp.22,36-38,536,541, 699-702.
- Reinherz EL, et al.: 1979. A monoclonal antibody with selective reactivity with functionally mature human thymocytes and all peripheral human T cells. J Immunol 123:1312-1317.
- Beverly PCL and Callard RE: 1981. Distinctive functional characteristics of human T lymphocytes defined by E rosetting or a monoclonal anti-T cells antibody. Eur J Immunol 11:329-334.
- National Committee for Clinical Laboratory Standards: Procedures for the collection of diagnostic blood specimen by venipuncture; Approved Standard. 1991. NCCLS Document H3-A3.
- Loken MR, et al.: 1990. Establishing optimal lymphocyte gates for immunophenotyping by flow cytometry. Cytometry 11:453-459.
- Nicholson JKA, et al.: 1996. Use of CD45 fluorescence and side scatter characteristics for gating lymphocytes when using the whole blood lysis procedure and flow cytometry. Cytometry 11:453-459.
- Gebel HM, et al.: 1989. Discordant expression of CD3 and T-cell receptor antigens on lymphocytes from patients treated with OKT3. Transplantation Proceedings 1:1745-1746.
- Tunnacliffe A, et al.: 1989. The majority of CD3 epitopes are conferred by the epsilon chain. In Leukocyte Typing IV, White Cell Differentiation Antigens. W Knapp et al., eds, Oxford: Oxford University Press. p. 295-296.
- Schroeder TI and Chatenoud L: Immunological monitoring during treatment with OKT3. Presented under the auspices of the American Society of Transplant Surgeons. Citation to be filled in when available from authors. Koepke JA and Landay AL: 1989. Precision and accuracy of absolute lymphocyte counts. Clin Immunol Immunopathol 52:19-27.
- Koepke JA and Landay AL: 1989. Precision and accuracy of absolute lymphocyte counts. Clin Immunol Immunopathol 52:19-27.

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