

**For Research Use Only. Not for use in diagnostic procedures.**

### SPECIFICITY

The CD34 glycoprotein is a transmembrane single chain molecule. Its molecular weight is about 110 kDa. The extracellular domain is heavily N- and O-glycosylated (1 – 4) and the cytoplasmic sequence reveals two sites for activated protein kinase C phosphorylation and one site for tyrosine phosphorylation (3). This antigen is the earliest known marker in human for hematopoietic progenitor cells (5, 6).

There are three classes of CD34 epitopes defined by differential sensitivity to enzymatic cleavage with glycoprotease from *Pasteurella haemolytica*, with chymopapain and with neuraminidase.

The Immu133 monoclonal antibody recognizes specifically a Class I epitope, neuraminidase sensitive, and glycoprotease sensitive which is regularly expressed on normal samples.

The CD34 molecule is expressed on virtually all hematopoietic precursor cells (7), including the multipotent stem cells (8). The CD34 molecule is the earliest marker in human for precursors of colony forming cells in the bone marrow (5, 6).

The CD34 glycoprotein is not restricted to hematopoietic progenitors (9) and has been detected on capillary endothelial cells (9 – 11), and on bone marrow stromal cells and their precursors (12).

The Immu133 monoclonal antibody has been assigned to the CD34 cluster of differentiation at the Vth International Workshop on Human Leukocyte Differentiation Antigens in Boston, USA, in 1993 (13).

### REAGENT

IOTest CD34 - PC5 Conjugated Antibody  
PN IM2649 – Liquid – 100 tests –  
10 µL / test.

<b>Clone</b>	Immu133
<b>Isotype</b>	IgG1 mouse
<b>Immunogen</b>	Cells from KG1a and TF1 cell lines
<b>Hybridoma</b>	P3-X63-Ag.8.653 x Balb/c spleen cells
<b>Source</b>	Ascites fluid
<b>Purification</b>	Ion exchange or affinity chromatography

**Conjugation** R-phycoerythrin-Cyanine 5 (PC5)

**Molar Ratio** PC5 / Ig : 0.7 – 1

**Fluorescence** Excites at 486–580 nm  
Emits at 660–680 nm

### REAGENT CONTENTS

This antibody is provided in phosphate-buffered saline pH 7.4, containing 0.1% sodium azide and 2 mg/mL bovine serum albumin.

### APPLICATION

Flow cytometry.

### STATEMENT OF WARNINGS

1. This reagent contains 0.1% sodium azide. Sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, wash excessively with water.
2. Specimens, samples and all material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
3. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
4. Do not use antibody beyond the expiration date on the label.
5. Do not expose reagents to strong light during storage or incubation.
6. Avoid microbial contamination of reagents or incorrect results might occur.
7. Use good laboratory practices when handling this reagent.

### STORAGE CONDITIONS AND STABILITY

This reagent is stable up to the expiration date when stored at 2 – 8°C in the dark. Do not freeze.

### REAGENT PREPARATION

No reconstitution is necessary. This monoclonal antibody may be used directly from the vial. Bring reagent to 18 – 25°C prior to use.

### PROCEDURE

This reagent is designed for flow cytometry. A wash is required to yield optimal results  
Assay volume: 10 µL per 5 x 10<sup>5</sup> cells in one test, or per 100 µL whole blood.

### DIRECT LABELLING OF WHOLE BLOOD OR BONE-MARROW WITH CONJUGATED MONOCLONAL ANTIBODY FOLLOWED BY RED BLOOD CELL LYSING METHOD.

Preliminary remarks:

In the following procedure, we recommend to use an ammonium chloride-based lysing solution with satisfactory results, provided that the sample is thoroughly washed with phosphate-buffered saline.

Procedure:

1. Pipet 100 µL of specimen into two tubes (control and test).
2. Add 30 µL of blocking mouse Ig at 300 µg/mL (optional).
3. Add 10 µL of CD34-PC5 to the test tube or 10 µL of conjugated isotypic control to the control tube.
4. Mix gently and incubate for 15 min. at room temperature (18-25 °C) in the dark.
5. Wash the preparation by adding 3 mL of PBS.
6. Centrifuge at 300 x g for 10 min. at room temperature.
7. Discard the supernatant and resuspend the pellet in 100 µL of PBS.
8. Lyse the stained specimen, according to the manufacturer's recommendations.
9. Wash the preparation by fulfilling the tubes with PBS.
10. Centrifuge at 300 x g for 10 min. at room temperature.
11. Discard the supernatant and resuspend the pellet in 1 mL of PBS + 0.5% formaldehyde.

The preparations are ready to be analyzed by flow cytometry within 2 hours.

(Keep the preparation at 2 to 8 °C for an analysis within 24 hours.)

### DIRECT LABELLING OF ISOLATED MONONUCLEAR CELLS OR CULTURED CELLS.

After isolation of mononuclear cells (MNC) by density gradient centrifugation:

Resuspend cells in PBS containing 0.5% BSA and 0.05% NaN<sub>3</sub>; count the viable cells (Trypan Blue test), and adjust the cell concentration to 5 x 10<sup>6</sup> to 1 x 10<sup>7</sup> cells/mL.  
Cell culture adjustment:

Resuspend the pellet in PBS containing 0.5% BSA and 0.05% NaN<sub>3</sub> and adjust cell concentration to 5 x 10<sup>6</sup> to 1 x 10<sup>7</sup> cells/mL.

Procedure for direct immunolabelling:


1. Pipet 100 µL of the suspension into two tubes (test and control).
2. Add 30 µL of blocking mouse Ig at 300 µg/mL (optional).
3. Add 10 µL of CD34-PC5 to the test tube, or 10 µL of conjugated isotypic control to the control tube.
4. Mix gently and incubate during 30 min. at 2-8 °C in the dark.
5. Wash the preparation by adding 3 mL of cold (2-8 °C) PBS containing 0.5% BSA and 0.05% NaN<sub>3</sub>.
6. Centrifuge at 300 x g for 5 min. at 2-8 °C.
7. Repeat steps number 5 and 6 two times.
8. Resuspend stained pellets in 1 mL of cold (2-8 °C) PBS. Analyze by flow cytometry within 2 hours.  
(Cells can be fixed in PBS containing 0.5% formaldehyde, and kept at 2-8 °C for analysis within 24 hours.)

## PN IM2649 – 100 tests – Liquid – 10 µL / test – Clone Immu133

### EXAMPLE DATA

The histogram below is a monoparametric representation (Count versus Fluorescence Intensity) of KG-1a cell line. Staining is with CD34-PC5 monoclonal antibody (see catalog for PN). The isotopic control (see catalog for PN) labeling is shown in light.

Acquisition is with a COULTER<sup>®</sup> EPICS<sup>®</sup> XL<sup>™</sup> flow cytometer. Analysis is with the XL<sup>™</sup> SYSTEM II<sup>™</sup> software.



### SELECTED RESEARCH REFERENCES

1. Watt, S.M., Karhi, K., Gatter, K., Furley, A.J.W., Katz, F.E., Healy, L.E., Altass, L.J., Bradley, N.J., Sutherland, D.R., Levinsky, R., Greaves, M.F., "Distribution and epitope analysis of the cell membrane glycoprotein (HPCA-1) associated with human hematopoietic progenitor cells", 1987, *Leukemia*, 1, 1, 417-426.
2. Sutherland, D.R., Watt, S.M., Dowder, G., Karhi, K., Baker, M.A., Greaves, M.F., Smart, J.E., "Structural and partial amino-acid sequence analysis of the human hematopoietic progenitor cell antigen CD34", 1988, *Leukemia*, 12, 2, 793-803.
3. Simmons, D.L., Satterthwaite, A.B., Tenen, D.G., Seed, B., "Molecular cloning of a cDNA encoding CD34: A sialomucin of human hematopoietic stem cells", 1992, *J. Immunol.*, 1, 148, 267-271.
4. Greaves, M.F., Brown, J., Molgaard, H.V., Spurr, N.K., Robertson, D., Delia, D., Sutherland, D.R., "Molecular features of CD34: a hemopoietic progenitor cell-associated molecule", 1992, *Leukemia*, 6, 31-36.
5. Andrews, R.G., Singer, J.W., Bernstein, I.D., "Human hematopoietic precursors in long-term culture, Single CD34+ cells that lack detectable T cell, B cells and myeloid cell antigens produce multiple colony-forming cells when cultured with marrow stromal cells", 1990, *J. Exp. Med.*, 172, 355-358.
6. Lansdorp, P.M., Sutherland, H.J., Eaves, C.J., "Selective expression of CD45 isoforms on functional subpopulations of CD34+ hemopoietic cells from human bone marrow", 1990, *J. Exp. Med.*, 172, 363-366.
7. Civin, C.I., Strauss, L.C., Brovall, c., Jackler, M.J., Schwartz, J.F., Shaper, J.H., "Antigenic analysis of hematopoiesis. III. A hematopoietic progenitor cell surface antigen defined by a monoclonal antibody raised against KG-1a cells", 1984, *J. Immunol.*, 1, 133, 157-165.
8. Terstappen, L.W.M.M., Huang, S., Safford, M., Lansdorp, P.M., Loken, M.R., "Sequential generation of hematopoietic colonies derived from single nonlineage-committed CD34+CD38- progenitor cells", 1991, *Blood*, 6, 77, 1218-1227.
9. Civin, C.I., Trischmann, T.M., Fackler, M.J., Bernstein, I.D., Bühring, H.J., Campos, L., Greaves, M.F., Kamoun, M., Katz, D.R., Lansdorp, P.M., Look, A.T., Seed, B., Sutherland, D.R., Tindle, R.W., Uchanska-Ziegler, B., "Summary of CD34 cluster workshop section", 1989, *Leucocyte Typing IV, White Cell Differentiation Antigens*. W. Knapp, et al., Eds., Oxford University Press, 818-825.
10. Delia, D., Lampugnani, M.G., Resnati, M., Dejana, E., Ajello, A., Fontanella, E., Soligo, D., Pierotti, M.A., Greaves, M.F., "CD34 expression is regulated reciprocally with adhesion molecules in vascular endothelial cells in vitro", 1993, *Blood*, 4, 81, 1001-1008.
11. Fina, L., Molgaard, H.V., Robertson, d., Bradley, N.J., Monaghan, P., Delia, D., Sutherland, D.R., Baker, M.A., Greaves, M.F., "Expression of the CD34 gene in vascular endothelial cells", 1990, *Blood*, 12, 75, 2417-2426.
12. Simmons, P.J., Torok-Storb, B., "CD34 expression by stromal precursors in normal human adult bone marrow", 1991, *Blood*, 78, 2848-2853.
13. Greaves, M.F., Titley, I., Colman, S.M., Bühring, H.-J., Campos, L., Castoldi, G.L., Garrido, F., Gaudernack, G., Girard, J.-P., Inglès-Esteve, J., Invernizzi, R., Knapp, W., Lansdorp, P.M., Lanza, F., Merle-Béral, H., Parravicini, C., Razak, K., Ruiz-Cabello, F., Springer, T.A., van der Schoot, C.E., Sutherland, D.R., "CD34 cluster Workshop report", 1995, *Leucocyte Typing V, White Cell Differentiation Antigens*. Schlossman, S.F., et al., Eds., Oxford University Press, 840-846.

### TRADEMARKS AND PATENTS

Beckman Coulter, the Beckman Coulter logo, COULTER, COULTER EPICS, IOTest, SYSTEM II, XL, software are registered trademarks of Beckman Coulter, Inc.

PC5 is licensed under patents 4,542,104 and 4,520,110.

Cy5 is licensed under patents 4,981,977 and 5,268,486.

### MANUFACTURED BY:

Immunotech SAS,  
a Beckman Coulter Company  
130, avenue de Lattre de Tassigny  
B.P. 177 - 13276 Marseille Cedex 9  
France

For additional information in the USA, call 800-526-7694.

Outside the USA, contact your local Beckman Coulter representative.

[www.beckmancoulter.com](http://www.beckmancoulter.com)

Printed in France.

Made in France.

©2008 Beckman Coulter, Inc.

All Rights Reserved