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1. Szollosi J, Damjanovich S, Matyus L, Application of fluorescence resonance energy transfer in the clinical laboratory: routine and research. *Cytometry* 1998; 34:159-179.
2. Didenko VV. DNA probes using fluorescence resonance energy transfer (FRET): designs and applications. *Biotechniques* 2001; 31(5): 1106-1121.
3. Chan FK, Siegel RM, Zacharias D, Swofford R, Holmes KL, Tsien RY et al. Fluorescence resonance energy transfer analysis of cell surface receptor interactions and signaling using spectral variants of the green fluorescent protein. *Cytometry* 2001; 44(4): 361-368.
4. Ashcroft RG, Lopez PA. Commercial high-speed machines open new opportunities in high throughput flow cytometry. *Journal of Immunological Methods* 2000; 243: 13-24.
5. Roncarati R, Sestan N, Scheinfeld MH, Berechid BE, Lopez PA, Meucci O et al. The  $\alpha$ -secretase-generated intracellular domain of  $\beta$ -amyloid precursor protein binds Numb and inhibits Notch signaling. *PNAS* 2002; 99(10): 7102-7107.
6. Scheinfeld MH, Roncarati R, Vito P, Lopez PA, Abdallah M, D'Adamio L. Jun NH2-terminal kinase (JNK) interacting protein 1 (JIP1) binds the cytoplasmic domain of the Alzheimer's  $\beta$ -amyloid precursor protein (APP). *J Biol Chem* 2002; 277(5): 3767-3775.

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