

**IOTest®**  
**CD14+16-FITC /**  
**CD85k-PE /**  
**CD123-PC5**

**REF A23416**

50 tests; 1 mL  
 20 µL / test



**IOTest**  
**Conjugated Antibodies**



ENGLISH	Specification of constituent 1	Specification of constituent 2	Specification of constituent 3
<b>Specificity</b>	CD14+16	CD85k (ILT3)	CD123
<b>Clone</b>	RMO52+3G8	ZM3.8	SSDCLY107D2
<b>Hybridoma</b>	SP2/0 x Balb/c	Ag8.653 x Balb/c	SP2/0 x Balb/c
<b>Immunogen</b>	Human monocytes+ neutrophils	Fusion protein of human ILT3	Plasmacytoid human Dendritic Cells
<b>Immunoglobulin</b>	IgG2a + IgG1	IgG1	IgG1
<b>Species</b>	Mouse	Mouse	Mouse
<b>Source</b>	Ascites	Ascites	Ascites
<b>Purification</b>	Protein A affinity chromatography	Protein A affinity chromatography	Protein A affinity chromatography
<b>Fluorochrome</b>	Fluorescein isothiocyanate (FITC)	Phycoerythrin (PE)	Phycoerythrin Cyanin 5.1 (PC5)
<b>λ excitation</b>	488 nm	488 nm	488 nm
<b>Emission peak</b>	525 nm	575 nm	670 nm
<b>Buffer</b>	PBS pH 7.2 plus 2 mg / mL BSA and 0.1% NaN <sub>3</sub>		

**USE**

This fluorochrome-conjugated antibody mixture is suitable for multiparametric analysis using flow cytometry. It permits the identification and enumeration of the plasmacytoid subset of circulating dendritic cells (PDC) in human biological samples.

**PRINCIPLE**

This test is based on the ability of specific monoclonal antibodies to bind to the antigenic determinants expressed by leucocytes. Specific staining of the leucocytes is performed by incubating the sample with the IOTest reagent. The red cells are then removed by lysis and the leucocytes, which are unaffected by this process, are analyzed by flow cytometry. The flow cytometer measures light diffusion and the fluorescence of cells. It makes possible the delimitation of the population of interest within the electronic window defined on a histogram (see examples in the Appendix for the description of the proposed gating strategy).

**EXAMPLES OF CLINICAL APPLICATIONS**

The simultaneous analysis of CD14, CD16, CD85k (ILT3) and CD123 antigens, helps the identification and the characterization of peripheral blood plasmacytoid dendritic cells (PDCs), based on the selection of CD123 expressing cells (basophils, PDCs and some monocytes), and the positive selection of PDCs + monocytes with ILT3. The discrimination of PDCs from monocyte subsets is then performed with CD14 + CD16. Recent studies point to a numerical decrease and sometimes even functional impairment of circulating plasmacytoid and/or myeloid dendritic cells subsets in various pathologies. Due to the high levels of IFNα they secrete, PDCs have been postulated to play a role in the pathogenesis of viral infections and Lupus (1). Studies have shown a decrease in the number of circulating PDCs in chronic Hepatitis B (2) or C patients (3) compared to non chronic and control individuals. Furthermore, in AIDS and cancer patients, the level of circulating DCs tends to decrease. Interestingly, a study in AIDS patients showed that, upon tri-therapy, the decrease of viral load correlated with an increase of circulating PDCs (4). In haematopoietic cancer patients, DCs counts may also be significantly reduced in lymphoid or myeloid leukaemia (5 - 6). A similar observation is found in certain solid cancers (7). Furthermore, circulating DCs are decreased in metastatic cancers as compared to patients with localized cancer (8).

**STORAGE AND STABILITY**

The conjugated liquid forms must be kept at between 2 and 8°C and protected from light, before and after the vial has been opened. Stability of closed vial: see expiry date on vial. Stability of open vial: the reagent is stable for 90 days.

**REAGENT CONTENTS**

Contact Beckman Coulter Customer Service to obtain the antibody concentration in the IOTest reagent.

**EVIDENCE OF DETERIORATION**

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address : immuno-techsup@beckmancoulter.com.

**PRECAUTIONS**

1. Do not use the reagent beyond the expiry date.
2. Do not freeze.
3. Let it come to room temperature (18 – 25°C) before use.
4. Minimize exposure to light.
5. Avoid microbial contamination of the reagents, or false results may occur.
6. Antibody solutions containing sodium azide (NaN<sub>3</sub>) should be handled with care. Do not ingest and avoid contact with the skin, mucosa and eyes. Moreover, in an acid medium, sodium azide can form the potentially dangerous hydrazoic acid. If it needs to be disposed of, it is recommended that the reagent be diluted in a large volume of water before pouring it into the drainage system so as to avoid the accumulation of sodium azide in metal pipes and to prevent the risk of explosion.
7. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
8. Never pipette by mouth and avoid all contact of the samples with the skin, mucosa and eyes.
9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.

**SAMPLES**

Venous blood or bone marrow samples must be taken using sterile tubes containing an EDTA salt as the anticoagulant. The use of other anticoagulants is not recommended. The samples should be kept at room temperature (18 – 25°C) and not shaken. The samples should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

Note: optimal results are obtained when the processed samples are analyzed within 6 hours after immunostaining.

**METHODOLOGY**

**NECESSARY MATERIAL NOT SUPPLIED**

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 10, 20, 100 and 500 µL.
- Plastic haemolysis tubes.
- Calibration beads. For example: Flow-Set™ Fluorospheres (Ref. 6607007).
- Red cell lysis reagent with washing stage after lysis. For example: VersaLyse™ (Ref. A09777).
- Leucocyte fixation reagent. For example : IOTest 3 Fixative Solution (Ref. A07800).
- Isotypic controls: IOTest reagents : IgG1-PC5 (Ref. A07798) and IgG1-FITC/IgG1-PE (Ref. A07794).
- Buffer (PBS: 0.01 M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

**PROCEDURE**

**NOTE:** The procedure below is valid for standard applications. Sample and/or VersaLyse volumes for certain Beckman Coulter applications may be different. If such is the case, follow the instructions on the application's technical leaflet. For each sample analyzed, in addition to the test tube, one control tube is required in which the cells are mixed in the presence of the isotypic controls.

1. Add 20 µL of specific IOTest conjugated antibodies to each test tube, and 20 µL of IgG1-FITC/IgG1-PE and 10 µL of IgG1-PC5 to the control tube.
2. Add 100 µL of the test sample to the 2 tubes. Vortex the tubes gently.
3. Incubate for 15 to 20 minutes at room temperature (18 – 25°C), protected from light.
4. Then perform lysis of the red cells, if necessary, by following the recommendations of the lysis reagent used. As an example, if you wish to use VersaLyse (Ref. A09777), refer to the leaflet and follow preferably the procedure called "with concomitant fixation", which consists of adding 1 mL of the "Fix-and-Lyse" mixture prepared extemporaneously. Vortex immediately for one second and incubate for 10 minutes at room temperature, protected from light.

If the sample does not contain red cells, add 2 mL of PBS.

5. Centrifuge for 5 minutes at 150 x g at room temperature.
6. Remove the supernatant by aspiration.
7. Resuspend the cell pellet using 3 mL of PBS.
8. Repeat step 5.
9. Remove the supernatant by aspiration and resuspend the cell pellet using:
  - 0.5 mL or 1 mL of PBS plus 0.1% of formaldehyde if the preparations are to be kept for more than 2 hours and less than **6 hours**. (A 0.1% formaldehyde PBS can be obtained by diluting 12.5 µL of the IOTest 3 Fixative Solution (Ref. A07800) at its 10X concentration in 1 mL of PBS).
  - 0.5 mL or 1 mL of PBS without formaldehyde, if the preparations are to be analyzed within 2 hours.

**NOTE:** In all cases, keep the preparations between 2 and 8°C and protected from light.

## PERFORMANCE

Performance data are obtained using the procedure described above on 24 hour-old blood samples previously collected on sterile tubes with EDTA salt as anticoagulant. Analysis is performed within 2 hours following immunostaining.

## SPECIFICITY

The CD14 molecule is found on cells of myelomonocytic lineage. It is strongly expressed on monocytes, macrophages, and weakly on neutrophils (9, 10). It is also present on pleural phagocytic cells and on reticular dendritic cells, on Langerhans cells, and histiocytes (11, 12). The monoclonal antibody (mAb) RMO52 does not react with T or B lymphocytes (9, 10).

MAb RMO52 was assigned to CD14 during the 6<sup>th</sup> HLDA Workshop on Human Leucocyte Differentiation Antigens, held in Kobe, Japan, in 1996 (WS Code: MA62, Section M) (13).

The CD16 antigen is the low-affinity receptor for IgG (FcγRIII). The CD16 antigen exists in two different forms encoded by two different genes: FcγRIIIA (or III-2) and FcγRIIIB (or III-1). One is a transmembrane form (FcγRIIIA) expressed on NK cells, monocytes and macrophages. The other is a glycosyl-phosphatidylinositol (GPI)-anchored form (FcγRIIIB) only expressed on neutrophils (14, 15).

MAb 3G8 binds to the FcγRIIIA as well as (strongly) to the FcγRIIIB.

MAb 3G8 was assigned to CD16 during the 4<sup>th</sup> HLDA Workshop of Vienna, Austria, in 1989 (WS Code: 409, Section NL) (16).

CD85k, alias ILT3, is a member of a family of genes located on human chromosome 19. The ILT gene-derived proteins are the Immunoglobulin (Ig)-Like Transcripts (ILT), also known as Leucocyte Ig-like Receptors (LIR) and Monocyte / Macrophage Ig-like Receptors (MIR) (17 – 19). ILT / LIR / MIR are preferentially expressed on monocytes, macrophages, dendritic cells (DCs) and granulocytes.

CD85k, otherwise known as ILT3, LIR-5, HM18, is a transmembrane protein of 60 kDa, with 2 extracellular Ig-SF domains, and 3 cytoplasmic ITIMs. The molecule is constitutively phosphorylated. It is selectively expressed by monocytes, macrophages and DCs (20).

MAb ZM3.8 is specific for human ILT3, it immunoprecipitates a protein of 60 kDa from monocytes, and is expressed in human blood samples by a subset of dendritic cells (20 - 22).

The CD123 antigen, or Interleukin-3 alpha receptor (IL-3R alpha), belongs to the cytokine receptor family (23, 24). It is constitutively expressed by committed hematopoietic stem / progenitor cells, by most of the myeloid lineage (CD13<sup>+</sup>, CD14<sup>+</sup>, CD33<sup>+</sup>, CD15<sup>low</sup>), and by some CD19<sup>+</sup> cells (25); but it is absent from CD3<sup>+</sup> cells (25). MAb SSDCLY107D2 reacts with CD123 molecule present on PDCs (23).

## LINEARITY

To test the linearity of staining for the CD123 and CD85k specificities of this reagent, a sample containing positive cells (PDC sorted from normal whole blood) and a negative cell line MOLTA (CD123<sup>-</sup> CD85k<sup>-</sup>) were mixed in different proportions and a constant final number of cells, so that the positive cells/negative cell line ratio of the mixture ranged from 0,02 to 0,15%. This corresponds to concentrations ranging from normal values to 10 fold less than normal values.

Aliquots were stained using the procedure described above and linear regression between the expected values and the observed values was calculated.

Specificity	Linear regression	Linearity (R <sup>2</sup> )
CD123 <sup>+</sup> CD85k <sup>+</sup>	Y = 1.0201 X - 0.0012	0.9953

## EXPECTED VALUES

Each laboratory must compile a list of reference values based upon a group of healthy donors from the local population. This must be done by taking age, sex and ethnic group into account, as well as any other potential regional differences.

In our laboratories, the whole blood samples of 30 healthy adults were used. The results obtained for the mean percentage of PDC amongst PBMC (peripheral blood mononuclear cells) are given in the table below:

PDC	n	Mean %	SD	CV (%)
CD14 <sup>-</sup> /CD16 <sup>-</sup> CD123 <sup>+</sup> CD85k <sup>+</sup>	30	0.495	0.203	40.99

## INTRA-LABORATORY REPRODUCIBILITY

On the same day and using the same cytometer, 12 measurements of the positivity of a sample containing positive cells (peripheral blood from the same donor) were carried out. The results obtained are summarized in the following table:

PDC	n	Mean %	SD	CV (%)
CD14 <sup>-</sup> /CD16 <sup>-</sup> CD123 <sup>+</sup> CD85k <sup>+</sup>	12	0.283	0.018	6.20

## LIMITATIONS OF THE TECHNIQUE

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence spillover have not been correctly compensated for and if the regions have not been carefully positioned.
2. It is preferable to use a RBC lysis technique with a washing step as this reagent has not been optimized for "no wash" lysis techniques.
3. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
4. The conjugated antibodies of this reagent are calibrated so as to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
5. In the case of a hyperleucocytosis, dilute the specimen in PBS so as to obtain a value of approximately 5 x 10<sup>9</sup> leucocytes/L.
6. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining.
7. Due to the tandem structure of the fluorochrome, PC5 also emits light at 575 nm. This secondary emission peak varies from lot-to-lot of PC5. Therefore, for multi-color analysis, the compensation matrix should be carefully checked when changing the lot of a PC5-conjugate.

## MISCELLANEOUS

See the Appendix for examples and references.

## TRADEMARKS

The Beckman Coulter logo, Beckman Coulter, COULTER, EPICS, Flow-Set, IOTest, System II, VersaLyse and XL are the registered trademarks of Beckman Coulter Inc.

## MANUFACTURED BY :

IMMUNOTECH SAS  
A Beckman Coulter Company  
130 avenue de Lattre de Tassigny  
B.P. 177 – 13276 Marseille Cedex 9  
France  
Customer Services: (33) 4 91 17 27 27

Cy5 is a trademark of Amersham Biosciences, a GE Healthcare company.

Cy5 is licensed under patents 4,981,977 and 5,268,486.

[www.beckmancoulter.com](http://www.beckmancoulter.com)

Printed in France.

Made in France.


© 2008 Beckman Coulter, Inc.  
All Rights Reserved.




## APPENDIX TO REF A23416

### EXAMPLES


The histograms below are biparametric representations of a lysed (see procedure described above) normal whole blood. The staining is performed with the following 3-color combination: (CD14+CD16)-FITC / CD85k (ILT3)-PE / CD123-PC5 ( Ref. A23416).



**Histogram 1:**  
SS versus FS, displaying all events. Region A is drawn around all white blood cells in order to eliminate debris. A Peripheral Blood Mononuclear Cell (PBMC) region is drawn around the lymphocytic and monocytic populations in order to be used for final calculation of PDC's percentage.




**Histogram 2 (gated on A):**  
FL4 (CD123-PC5) versus SS, displaying all gated events. Region B is drawn in order to exclude all CD123-negative cells.



**Histogram 3 (gated on A AND B):**

FL1 (CD14+CD16)-FITC versus FL2 (ILT3-PE). Region C is drawn around all CD(14+16)<sup>neg</sup> / ILT3<sup>+</sup>. As a visual control of the gating, the eosinophils (Eo) should be visible in the bottom right corner of the region C (see "Specificity").



**Histogram 4 (gated on A AND B AND C):**

FL4 (CD123-PC5) versus FL2 (ILT3-PE). Region D is drawn around CD123<sup>bright</sup> / ILT3<sup>bright</sup> events. (the percentage of PDCs amongst PBMC is calculated by dividing the number of events in gate D by the number of events in gate PBMC as defined on histogram 1). Acquisition is with a COULTER® EPICS® XL™ flow cytometer, using System II™ acquisition software. Analysis is with Beckman Coulter CXP Cytometer software.

## REFERENCES

1. Palucka, A.K., Banchereau, J., Blanco, P., Pascual, V., "The interplay of dendritic cell subsets in systemic lupus erythematosus", 2002, *Immunol. Cell. Biol.*, 80, 484-488.
2. Duan, X.Z., Wang, M., Li, H.W., Zhuang, H., Xu, D., Wang, F.S., "Decreased Frequency and Function of Circulating Plasmacytoid Dendritic Cells (PDC) in Hepatitis B Virus Infected Humans", *J Clin Immunol.* 2004 Nov, 24, 6, 637-646.
3. Tsubouchi, E., Akbar, S.M., Horiike, N., Onji, M., "Infection and dysfunction of circulating blood dendritic cells and their subsets in chronic hepatitis C virus infection" *J. Gastroenterol.*, 2004, 8, 754-762.
4. Pacanowski, J., Kahi, S., Baillet, M., Lebon, P., Deveau, C., Goujard, C., Meyer, L., Oksenhender, E., Sinet, M., Hosmalin, A., "Reduced blood CD123+ (lymphoid) and CD11c+ (myeloid) dendritic cell numbers in primary HIV-1 infection", 2001, *Blood*, 98, 3016-3021.
5. Mami, N., B., Mohty, M., Chambost, H., Gaugler B., Olive, D., "Blood dendritic cells in patients with acute lymphoblastic leukaemia", 2004, *British Journal of Haematology*, 126, 1, 77 - 80.
6. Mohty, M., Isnardon, D., Vey, N., Brière, F., Blaise, D., Olive D., Gaugler, B., "Low blood dendritic cells in chronic myeloid leukaemia patients correlates with loss of CD34<sup>+</sup>/CD38<sup>+</sup> primitive haematopoietic progenitors", 2002, *British Journal of Haematology*, 119, 1, 115 - 118.
7. Yanagimoto, H., Takai, S., Sato, S., Toyokawa, H., Takahashi, K., Terakawa, N., Kwon, A.H., Kamiyama, Y., "Impaired function of circulating dendritic cells in patients with pancreatic cancer", *Clin Immunol.* 2005, 114, 1, 52-60.
8. Lissoni, P., Vigore, L., Ferranti, R., Bukovec, R., Meregalli, S., Mandala, M., Barni, S., Tancini, G., Fumagalli, L., Giani, L., "Circulating dendritic cells in early and advanced cancer patients: diminished percent in the metastatic disease", 1999, *J. Biol. Regul. Homeost. Agents*, 13, 216-219.
9. Todd III, R.F., Nadler, L.M., Schlossman, S.F., "Antigens on human monocytes identified by monoclonal antibodies", 1981, *J. Immunol.*, 4, 126, 1435-1442.
10. Todd, R.F., van Agthoven, A., Schlossmann, S.F., Terhorst, C., "Structural analysis of differentiation antigens Mo1 and Mo2 on human monocytes", 1982, *Hybridoma*, 3, 1, 329-337.
11. Peters, J.H., Ruppert, J., Gieseler, R.K.H., Najjar, H.M., Xu, H., "Differentiation of human monocytes into CD14 negative accessory cells: do dendritic cells derive from the monocytic lineage?", 1991, *Pathobiology*, 59, 122-126.
12. Ziegler-Heitbrock, H.W.L., Ulevitch, R.J., "CD14: Cell surface receptor and differentiation markers", 1993, *Immunol. Today*, 3, 14, 121-125.
13. Goyert, S.M., Cohen, L., Gangloff, S.C., Ashmun, R., Haeflner-Cavaillon, N., "CD14 Workshop panel report", 1997, *Leucocyte Typing VI, White Cell Differentiation Antigens*. Kishimoto, T., et al, Eds., Garland Publishing, Inc., 963-965.
14. Ravetch, J.V., Perussia, B., "Alternative membrane forms of FcγRIII (CD16) on human natural killer cells and neutrophils", 1989, *J. Exp. Med.*, 170, 481-497.
15. Huizinga, T.W.J., Roos, D., Von dem Borne, A.E.G. Kr., "Neutrophil FcγRIII receptors: a two-way bridge in the immune system", 1990, *Blood*, 75, 1211-1214.
16. Schmidt, R.E., Perussia, B., "Cluster report: CD16", 1989, *Leucocyte Typing IV, White Cell Differentiation Antigens*. W. Knapp, et al., Eds., Oxford University Press, 574-578.
17. Samaridis, J., and Colonna, M., "Cloning of novel immunoglobulin superfamily receptors expressed on human myeloid and lymphoid cells: structural evidence for new stimulatory and inhibitory pathways", 1997, *Eur. J. Immunol.*, 3, 27, 660-665.
18. André, P., Biassoni, R., Colonna, M., Cosman, D., Lanier, L.L., Long, O., Lopez-Botet, M., Moretta, A., Moretta, L., Parham, P., Trowsdale, J., Vivier, E., Wagtmann, N., and Wilson, M.J., "New nomenclature for MHC receptors", 2001, *Nature Immunol.*, 2, 8, 661.
19. Marsh, S.G.E., Parham, P., Dupont, B., Geraghty, D.E., Trowsdale, J., Middleton, D., Vilches, C., Carrington, M., Witt, C., Guethlein, L.A., Shilling, H., Garcia, C.A., Hsu, K.C., and Wain, H., "Killer-cell immunoglobulin-like receptor (KIR) nomenclature report, 2002", 2003, *Immunogenetics*, 55, 220-226.
20. Cella, M., Döhning, C., Samaridis, J., Dessing, M., Brockhaus, M., Lanzavecchia, A., Colonna, M., "A novel inhibitory receptor (ILT3) expressed on monocytes, macrophages, and dendritic cells involved in antigen processing", 1997, *J. Exp. Med.*, 185, 1743-1751.
21. Mohty, M., Jarrossay, D., Lafage-Pochitaloff, M., Zandotti, C., Brière, F., de Lamballeri, X.N., Isnardon, D., Sainty, D., Olive, D., Gaugler, B., "Circulating blood dendritic cells from myeloid leukemia patients display quantitative and cytogenetic abnormalities as well as functional impairment", 2001, *Blood*, 13, 98, 3750-3756.
22. Chaperot, L., Bendriss, N., Manches, O., Gressin, R., Maynadie, M., Trimoreau, F., Orfeuvre, H., Corront, B., Feuillard, J., Sotto, J.J., Bensa, J.C., Brière, F., Plumas, J., Jacob, M.C., "Identification of a leukemic counterpart of the plasmacytoid dendritic cells", 2001, *Blood*, 10, 97, 3210-3217.
23. Miyajima, A., "CDw123 (interleukin 3 receptor alpha chain) Workshop Panel Report", 1997, *Leucocyte Typing VI., White Cell Differentiation Antigens*. Kishimoto, T., et al, Eds., Garland Publishing, Inc., 854-855.
24. Moretti, S., Lanza, F., Dabusti, M., Tieghi, A., Campioni, D., Dominici, M. and Castoldi, G.L., "CD123 (Interleukin 3 receptor alpha chain)", 2001, *J. Regul. Homeost. Agents.*, 15, 98-100.
25. Sato, N., Caux, C., Kitamura, T., Watanabe, Y., Arai, K-I, Banchereau, J., and Miyajima, A., "Expression and factor-dependent modulation of the interleukin-3 receptor subunits on human hematopoietic cells", 1993, *Blood*, 3, 82, 752-761.