

# IOTest CD56-APC-Alexa Fluor 700

PN B10822 – 0.5 mL – Liquid – 10 µL/test – Clone N901 (NKH-1)

## Analyte Specific Reagent.

Analytical and performance characteristics are not established.

### SPECIFICITY

The molecular weight of the CD56 heavily glycosylated protein (NKH1 molecule) is 200 – 220 kDa (1, 2). It is expressed on a subpopulation of peripheral blood lymphocytes that demonstrate non-major histocompatibility complex restricted cytotoxicity (1, 3).

The N901 (NKH-1) monoclonal antibody reacts with the majority of NK cells (1, 2). It also reacts with a minor subpopulation of CD3 positive T cells that mediates reduced cytotoxic activity (3). This antibody does not react with monocytes, granulocytes, erythrocytes or B lymphocytes.

The N901 (NKH-1) antibody has been assigned to the CD56 cluster of differentiation during the 4th International Workshop on Human Leucocyte Differentiation Antigens in Vienna, Austria, in 1989 (4).

### REAGENT

IOTest CD56-APC-Alexa Fluor 700

Conjugated antibody

PN B10822 - 0.5 mL - Liquid - 10 µL/test

<b>Clone</b>	N901 (NKH-1)
<b>Isotype</b>	IgG1, Mouse
<b>Immunogen</b>	Chronic myeloid leukemia (hCML) cells
<b>Hybridoma Source</b>	NS1 x balb/c Ascites fluid or supernatant of in vitro cultured hybridoma cells.
<b>Purification</b>	Affinity chromatography
<b>Conjugation</b>	Allophycocyanin-Alexa Fluor 700 (APC-Alexa Fluor 700)
<b>Molar Ratio</b>	APC-Alexa Fluor 700 / Ig : 0.5 - 1.5
<b>Fluorescence</b>	Excites at 633/638 nm Emits at 720 nm

### REAGENT CONTENTS

This antibody is provided in phosphate-buffered saline, containing 0.1% sodium azide and 2 mg/mL bovine serum albumin.

### STATEMENTS OF WARNING

1. This reagent contains 0.1% sodium azide. Sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, wash excessively with water.

2. Specimens, samples and all material coming in contact with them should be considered potentially infectious and disposed of with proper precautions.
3. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
4. Do not use antibody beyond the expiration date on the label.
5. Do not expose reagents to strong light during storage or incubation.
6. Avoid microbial contamination of reagents or incorrect results might occur.
7. Use good laboratory practices when handling this reagent.

### STORAGE AND HANDLING CONDITIONS AND STABILITY

This reagent is stable up to the expiration date when stored at 2 – 8°C. Do not freeze.

No reconstitution is necessary. This monoclonal antibody may be used directly from the vial. Bring reagent to 18 – 25°C prior to use.

### PRECAUTIONS

Due to the tandem structure of the fluorochrome, APC-Alexa Fluor 700 also emits light at 660 nm. This secondary emission peak varies from lot-to-lot of APC-Alexa Fluor 700. Therefore, for multi-color analysis, the compensation matrix should be carefully checked when changing the lot of a APC-Alexa Fluor 700-conjugate.

Weak non-specific binding on a lymphocyte subpopulation may occur on some donors with APC-Alexa Fluor 700 conjugates.

### SELECTED RESEARCH REFERENCES

1. Griffin, J.D., Hercend, T., Beveridge, R., Schlossman, S.F., "Characterization of an antigen expressed by human natural killer cells", 1983, J. Immunol., 130, 2947-2951.
2. Hercend, T., Griffin, J.D., Bensussan, A., Schmidt, R.E., Edson, M.A., Brennan, A., Murray, C., Daley, J.F., Schlossman, S.F., Ritz, J., "Generation of monoclonal antibodies to a human natural killer clone: Characterization of two natural killer-associated antigens, NKH1A and NKH2, expressed on subsets of large granular lymphocytes", 1985, J. Clin. Invest., 75, 932-943.
3. Lanier, L.L., Le, A.M., Civin, C.I., Loken, M.R., Phillips, J.H., "The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes", 1986, J. Immunol., 136, 4480-4486.

4. Schubert, J., Lanier, L.L., Schmidt, R.E., "Cluster report: CD56", 1989, Leucocyte Typing IV, White Cell Differentiation Antigens. W. Knapp, et al., Eds., Oxford University Press, 699-702.

### TRADEMARKS

Beckman Coulter logo and IOTest, are trademarks of Beckman Coulter; Beckman Coulter logo, IOTest are registered in the USPTO and SIPO.

Alexa Fluor is a trademark of Molecular Probes, Inc.

### MANUFACTURED BY :

IMMUNOTECH SAS  
a Beckman Coulter Company  
130, avenue de Lattre de Tassigny  
B.P. 177 - 13276 Marseille Cedex 9  
France

For additional information in the USA, call 800-526-7694.

Outside the USA, contact your local Beckman Coulter representative.

[www.beckmancoulter.com](http://www.beckmancoulter.com)

Printed in France.

Made in France.

©2011 Beckman Coulter, Inc.