

For Research Use Only. Not For Use In Diagnostic Procedures.

SPECIFICITY

CD1a is a member of the lipid antigen-presenting molecules family (1). Like the other CD1 proteins (group 1: CD1a, CD1b, CD1c and group 2: CD1d, CD1e), CD1a is a transmembrane protein of 43 – 49 kDa, non-covalently associated with β 2-microglobulin (2). The structure and functions of CD1 molecules are reviewed in Ref. 3.

As other group 1 CD1 antigens, CD1a is found on professional antigen-presenting cells (APC) (3) such as dendritic cells and peripheral node mantle cells (3). CD1a is also strongly expressed by cortical thymocytes (4).

The BL6 monoclonal antibody has originally been described as specific for thymocytes (5) and further assigned to the CD1a cluster of differentiation at the 3rd International Workshop on Human Leucocyte Differentiation Antigens in Oxford, England, in 1989 (WS Code: 427) (6, 7).

REAGENT

IOTest CD1a-APC Conjugated Antibody
PN IM3645 - 100 tests – Liquid 10 – µL/test

Clone	BL6
Isotype	IgG1 Mouse
Immunogen	Human thymocytes
Hybridoma	SP2/0-Ag14 x Balb/c
Source	Ascites fluid
Purification	Affinity chromatography on protein A
Conjugation	Allophycocyanin (APC)
Molar Ratio	APC / Ig.0.5 – 1.5
Fluorescence	Excites at 633 – 635 nm Emits at 660 nm

APC conjugates are recommended for use only on flow cytometers equipped with an exciting source of 633 nm (He-Ne laser) or 635 nm (Red diode laser).

REAGENT CONTENTS

This antibody is provided in phosphate-buffered saline pH 7.4, containing 0.1% sodium azide and 2 mg/mL bovine serum albumin.

APPLICATION

Flow Cytometry

STATEMENT OF WARNINGS

1. This reagent contains 0.1% sodium azide. Sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin

or eye contact occurs, wash excessively with water.

2. Specimens, samples and all material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
3. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
4. Do not use antibody beyond the expiration date on the label.
5. Do not expose reagents to strong light during storage or incubation.
6. Avoid microbial contamination of reagents or incorrect results might occur.
7. Use good laboratory practises when handling this reagent.

STORAGE CONDITIONS AND STABILITY

This reagent is stable up to the expiration date when stored at 2 – 8°C in the dark. Do not freeze.

REAGENT PREPARATION

No reconstitution is necessary. This monoclonal antibody may be used directly from the vial. Bring reagent to 18 – 25°C prior to use.

PROCEDURE

This reagent is designed for Flow Cytometry. Assay volume: 10 µL per 5 x 10⁵ cells in one test, or per 100 µL whole blood. A wash is required to yield optimal results.

SELECTED RESEARCH REFERENCES

1. Porcelli, S.A., Segelke, B.W., Sugita, M., Wilson, I.A, Brenner, M.B., "The CD1 family of lipid antigen-presenting molecules", 1998, *Immunol. Today*, 8, 19, 362-368.
2. Martin, L.H., Calabi, F., Lefebvre, F.A., Bilsland, C.A.G., Milstein, C., "Structure and expression of the human thymocyte antigens CD1a, CD1b, and CD1c", 1987, *Proc. Natl. Acad. Sci. USA*, 84, 9189-9193.
3. Melian, A., Beckman, E.M., Porcelli, S.A., Brenner, M.B., "Antigen presentation by CD1 and MHC-encoded class I-like molecules", 1996, *Curr. Opin. Immunol.*, 8, 82-88.
4. Ito, T., Inaba, M., Inaba, K., Toki, J., Sogo, S., Iguchi, T., Adachi, Y., Yamaguchi, K., Amakawa, R., Valladeau, J., Saeland, S., Fukuhara, S., Ikehara, S., "A CD1a⁺/CD11c⁺ subset of human blood dendritic cells is a direct precursor of Langerhans cells", 1999, *J. Immunol.*, 163, 1409-1419.
5. Grassi, F., Dezutter-Dambuyant, C., McIlroy, D., Jacquet, C., Yoneda, K., Imamura, S., Boumsell, L., Schmitt, D.,

- Autran, B., Debré, P., Hosmalin, A, "Monocyte-derived dendritic cells have a phenotype comparable to that dermal dendritic cells and display ultrastructural granules distinct from Birbeck granules", 1998, *J. Leucocyte Biol.*, 64, 484-493
6. Boumsell, L., "Cluster report: CD1", 1989, *Leucocyte Typing IV, White Cell Differentiation Antigens*, W. Knapp, et al., Eds., Oxford University Press, 251-254.
 7. Kahn-Perles, B., Wietzerbin, J., Caillol D.H., Lemonnier, F., "Delineation of three subsets of class I human T antigens (HTA) on Molt-4 cells, serologic and regulatory relationship to HLA class I antigens", 1985, *J. Immunol.*, 3, 134, 1759-1765.
 8. Boumsell, L., Knowles, R., "Summary of the T1 workshop", 1987, *Leucocyte Typing III, White Cell Differentiation Antigens*, McMichael A.J., et al., Eds., Oxford University Press, 71-72.
 9. McMichael, A.J., Gotch, F.M., "T-cell antigens: New and previously defined clusters", 1987, *Leucocyte Typing III, White Cell Differentiation Antigens*, McMichael A.J., et al., Eds., Oxford University Press, 31-62.
 10. Chaperot, L., Bendriss, N., Manches, O., Gressin, R., Maynadie, M., Trimoreau, F., Orfeuvre, H., Corront, B., Feuillard, J., Sotto, J.J., Bensa, J.C., Brière, F., Plumas, J., Jacob, M.C., "Identification of a leukemic counterpart of the plasmacytoid cells", 2001, *Blood.*, 97, 3210-3217.

TRADEMARKS AND PATENTS

Beckman Coulter, The Beckman Coulter logo, COULTER and IOTest are trademarks of Beckman Coulter, Inc.

MANUFACTURED BY:

Immunotech SAS,
A Beckman Coulter Company
130, avenue de Lattre de Tassigny, B.P. 177
13276 Marseille Cedex 9, France

For additional information in the USA, call 800-526-7694.

Outside the USA, contact your local Beckman Coulter representative.

www.beckmancoulter.com

Printed In France
Made in France

©2008 Beckman Coulter, Inc.
All Rights Reserved